



# Ready-to-use mRNA, delivered with speed and certainty

## High-quality mRNA means more than meeting purity thresholds

Hongene is the home of nucleic acids manufacturing — your specialist CDMO partner empowering the production of next-generation RNA therapeutics with integrated raw materials supply, end-to-end expertise, exceptional speed and a commitment to environmental sustainability.

Our off-the-shelf mRNA products offer a fast and reliable way to advance RNA research workflows. Following DNA template synthesis then IVT with co-transcriptional capping, each construct is purified and isolated with exceptional integrity, capping efficiency, and <0.01% dsRNA content, ensuring consistent and reproducible performance across studies and assays.

Choose from a growing portfolio of linear mRNA, saRNA, and circular RNA, or partner with our team to scale seamlessly to GMP-grade, LNP-formulated, or custom constructs enabled by our vertically integrated CDMO services.

Explore Hongene's full mRNA portfolio and services at

[hongene.com/mRNA](https://hongene.com/mRNA)

Available immediately from our California warehouse.

## From discovery to delivery

Hongene accelerates every stage of RNA innovation.

- 1. Comprehensive Raw Material supply**
- 2. Off-the-shelf and custom RNA (Research)**
- 3. Analytical & Process Development**
- 4. GMP CDMO services**


NTPs, cap analogs, enzymes, DNA

Method development, characterization, scale-up support

DNA synthesis, cationic lipids, IVT

LNP formulation

Fill & finish

 29520 Kohoutek Way,  
Union City, CA 94587

# Off-the-shelf mRNA portfolio

All products available for immediate shipment.



## Linear mRNA & saRNA

Catalog #	Product	Standard Aliquot
ON-242	EPO (N1-Me-pUTP)	1 mg
ON-244	Firefly Luciferase (Fluc) mRNA	1 mg
ON-245	Cas9 mRNA (N1-Me-pUTP)	1 mg
ON-278	EGFP mRNA (N1-Me-pUTP)	1 mg
ON-279	Fluc mRNA (N1-Me-pUTP)	1 mg
ON-297	COVID-19 Spike Protein Full Length-mRNA	1 mg
ON-301	saRNA EGFP	1 mg
ON-302	saRNA Fluc	1 mg
ON-480	Cre mRNA (N1-Me-pUTP)	1 mg
ON-481	mCherry mRNA (N1-Me-pUTP)	1 mg
ON-516	EGFP mRNA (5-OMe U)	1 mg
ON-540	Fluc mRNA (5-OMe U)	1 mg
ON-847	NanoLuciferase mRNA (N1-Me-Pseudo U)	1 mg
ON-848	OVA mRNA (N1-Me-Pseudo U)	1 mg
ON-972	Cy5-UTP labeled Fluc mRNA (N1-Me-pUTP)	1 mg

## Circular RNA

Cat. No.	Product	Standard Aliquot
ON-505	Fluc circRNA	0.5 mg
ON-506	EGFP circRNA	0.5 mg

Comprehensive raw  
materials supply


Off-the-shelf and custom  
RNA products

GMP CDMO services  
CDMO services (API and LNP  
formulation & fill/finish)

Don't see what you need? Please inquire about our custom synthesis services

[hongene.com/mRNA](https://hongene.com/mRNA)

[product.info@hongene.com](mailto:product.info@hongene.com)

 29520 Kohoutek Way,  
Union City, CA 94587

 **Hongene  
Biotech**